



**TECHNIMOUNT**  
EMS®

# WALL MOUNT PRO - FL

## USER MANUAL



SAFETY AND FLEXIBILITY  
WHERE IT MATTERS MOST

## Copyrights

Copyright 2025 Technimount E.M.S. Holding Inc. is an affiliated company of Technologies CGC Inc. All Rights Reserved. No part of this publication may be reproduced, transmitted, transcribed, backed-up on a retrieval system, translated in any language, modified in any form or by any means, without the written permission of Technimount or its affiliate companies.

## Trademarks

Product and medical device manufacturer names are Trademarks™ or Registered Trademarks® of their respective holders.

- Technimount EMS® and Bracket Pro Serie® are a registered trademarks of Technologies CGC Inc.
- Oxivir® is a registered trademark of Diversey.
- Lavo® is a registered trademark of KIK Holdco Company, Inc.

For any issues with your Technimount product, its components, or for any technical questions during the installation, operation, or maintenance, please contact Technical Support at [techsupport@technimount.com](mailto:techsupport@technimount.com).

## Contact Information

### Technimount E.M.S. Holding Inc.

C/O: Regulatory Affairs

3770 Jean Marchand Street, Suite 100-C

Quebec (Quebec) G2C 1Y6

Canada

[customerservice@technimount.com](mailto:customerservice@technimount.com)

[techsupport@technimount.com](mailto:techsupport@technimount.com)

[legal@technimount.com](mailto:legal@technimount.com)

[www.technimount.com](http://www.technimount.com)

T + 1 581 700-6613

TF + 1 888 639-2758 (North America)

## Table of Contents

<b>Copyrights</b> .....	<b>ii</b>
<b>Trademarks</b> .....	<b>ii</b>
<b>Contact Information</b> .....	<b>ii</b>
<b>1. General Mentions and Considerations</b> .....	<b>5</b>
..... 1.1. Intended Use .....	5
..... 1.2. User Competency .....	5
..... 1.3. Warranties .....	6
..... 1.3.1. Warranty Policy .....	6
..... 1.3.2. Limited Warranty .....	6
..... 1.3.3. International Warranty Clause .....	6
..... 1.3.4. User Liability .....	7
..... 1.4. Claims .....	7
..... 1.4.1. Damaged or Defective Merchandise .....	7
..... 1.4.2. Return Policy .....	7
..... 1.4.3. Return of Material Authorization (RMA) .....	8
..... 1.4.4. Claim Process .....	8
<b>2. General Safety Guidelines</b> .....	<b>9</b>
..... 2.1. Symbols and Definitions .....	9
..... 2.2. Labels .....	10
..... 2.3. Safety Measures .....	11
<b>3. Technical Specifications</b> .....	<b>12</b>
<b>4. Wall Mount Pro - FL Orientation Diagrams</b> .....	<b>13</b>
<b>5. Wall Mount Pro - FL Illustrated Parts</b> .....	<b>14</b>
<b>6. Operate the Wall Mount Pro - FL</b> .....	<b>16</b>
..... 6.1. Install a Mounting System on the Wall Mount Pro - FL .....	16
..... 6.2. Remove a Mounting System from the Wall Mount Pro - FL .....	18
<b>Annex I EMS and clinical personnel Skills Assessment</b> .....	<b>19</b>
<b>Annex II Unpack the Wall Mount Pro - FL</b> .....	<b>21</b>
<b>Annex III Install the Wall Mount Pro - FL</b> .....	<b>23</b>
..... Required Installation Time .....	23
..... Required Tools .....	23
..... Required Parts .....	24

**Annex IV Maintenance..... 27**

..... Maintenance Frequency.....28

..... Required Tools .....28

..... Tested Cleaning Solutions.....28

..... Maintenance Plan.....29

..... Illustrated Inspection Points.....31

**Annex V Replacement Parts/Kits ..... 33**

## 1. General Mentions and Considerations

The Wall Mount Pro - FL user manual includes detailed product information, standards and guidelines to assist the administrator/manager/supervisor and biomedical technician (or equivalent) with the unpacking, assembling (when indicated) and maintenance of the Technimount product. It also includes specific user-related information to assist trained EMS and clinical personnel with effectively operating the mounting system.

Please read the user manual thoroughly to fully assess, comprehend, then relay its content to EMS and clinical personnel during training, to warn them of any potential danger of its abuse, how to safely use the product and provide a safe environment for patients as well as themselves. Your existing protocols should be updated to include the Technimount product(s) standards, guidelines, requirements and safety recommendations included within this documentation. The user manual should remain available to users when needed and relayed if the product is subsequently sold.

**NOTE :** Technimount continually seeks advancements in product design and quality. While the user manual contains the most updated product information available at the time of printing, it may contain minor differences from the current version, including image references. For more information, please contact Technical Support at [techsupport@technimount.com](mailto:techsupport@technimount.com).

**NOTE :** Technimount reserves the right to change part numbers and products without notice. Please contact Customer Service at [customerservice@technimount.com](mailto:customerservice@technimount.com) to ensure product options and availability.

### 1.1. Intended Use

The Wall Mount Pro - FL is designed to aid trained EMS and clinical personnel to hold and secure on vertical surfaces in airplanes and helicopters, medical devices and accessories that are essential to patient care during air EMS transport only.

### 1.2. User Competency

To safely operate the mounting system, EMS and clinical personnel must have the required skill level. Training should be given to EMS and clinical personnel, taking in account the skill level that is necessary to comply with their function and level of interaction with the Wall Mount Pro - FL:

- **Proficient (trained EMS and clinical personnel):** Has received the required training, is sufficiently knowledgeable to safely operate the product and have passed the skills assessment (refer to « Annex I EMS and clinical personnel Skills Assessment » on page 19).

**NOTE :** Any member of the EMS and clinical personnel who has not received the required training and lacks the knowledge needed to safely operate the mounting system must not use the product.

- **Expert (administrator/manager/supervisor):** Has in-depth knowledge and product comprehension, and is familiar with standards and guidelines. Skilled to train EMS and clinical personnel on how to safely use the product.

- **Advanced (biomedical technician or equivalent):** Has extensive mechanical experience. Skilled to perform the unpacking, assembly, safety checks and condition-based maintenance procedures as detailed in « Annex IV Maintenance » on page 27, basic troubleshooting, upgrade procedures and replacement procedures.

### 1.3. Warranties

#### 1.3.1. Warranty Policy

This statement constitutes Technimount's entire warranty policy with regards to Technimount products. Technimount makes no other warranty or representation, neither expressed nor implied, except as stated herein. There is no warranty of merchantability or warranty of fitness for any particular purpose. Under no circumstances will Technimount be held liable hereunder for incidental or consequential damages, arising from or in any manner, related to sales or use of any such product.

Technimount E.M.S. Holding Inc. guarantees to the original "Purchaser" of the "Product" with which this "Limited warranty" is included, that the product will be free from "Defects" in workmanship and materials under normal use for a "Warranty period" of one (1) year from the product purchase date by the purchaser. During the warranty period, the product will be repaired or replaced according to the "Limited warranty" without charge to the purchaser for parts or labor. The parts and product may be repaired or replaced with new or refurbished parts or products. Herein this Limited Warranty, "Refurbished" means parts and products which have been returned to the factory, specifically. If the product is repaired or replaced within the warranty period, the greater of the remaining warranty period will apply, or three (3) months from the date of repair or replacement. If the product is repaired or replaced after the warranty period has expired, the warranty period for the repair or replacement will expire three (3) months after the repair or replacement date.

#### 1.3.2. Limited Warranty

Technimount products are intended to retain medical devices in place in the case of a single emergency landing. Technimount products must not be reused if involved in a single emergency landing and must thereafter be replaced. If the end user uses a Technimount product following a single emergency landing, it is at the end user's own risk and Technimount will not be held liable.

The limited warranty does not apply to normal wear that could result from normal use. It does not apply when the product or any of its components have been disassembled or repaired by someone not authorized by "Technimount". It does not cover repair or the replacement of any product or part thereof damaged by neglect, misuse, moisture, liquids, exposure to heat, accidents, abuse, and non-compliance with the instructions for installation and use provided with the product. It does not cover physical damage to the surface of the product. The decision to repair, replace or refuse the coverage is final and at the sole discretion of Technimount, without any compensation or obligation from Technimount. The product defined as a "mounting system" is specifically designed to fill this requirement. Any other use will void the warranty and Technimount shall not be held liable on any claim if the product has been modified or adapted for use.

#### 1.3.3. International Warranty Clause

This warranty abides by the Canadian domestic policy. Warranty outside Canada may vary by country. Please contact Customer Service at [customerservice@technimount.com](mailto:customerservice@technimount.com) for more information.

#### 1.3.4. User Liability

The purchaser and administrator are responsible to validate regulations and standards for safety in their region, to comply with applicable safety regulations. Technimount is not responsible to inform the purchaser or the administrator of any applicable legislation for safety in their area.

The administrator is responsible for providing proper training to any personnel who will install, operate and perform maintenance on Technimount products.

#### 1.4. Claims

##### 1.4.1. Damaged or Defective Merchandise

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of the reception date. **Do not** accept damaged merchandise unless such damage is noted on the delivery receipt at the time of reception. Upon prompt notification, Technimount will file a freight claim with the appropriate carrier for damages incurred. Claims are limited in amount to the actual replacement cost. If the claim has not been received by Technimount within the fifteen (15) day period following the date of delivery, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for the full payment of the original invoice.

Claims for any short or broken merchandise must be made within thirty (30) days of invoicing. For details, refer to the claim process or contact Customer Service at [customerservice@technimount.com](mailto:customerservice@technimount.com).

##### 1.4.2. Return Policy

Technimount products may be returned up to sixty (60) days from the reception date, if:

- The received product does not match what was originally ordered.
- The product does not meet the Technimount technical sheet specifications.
- The product is not compatible with the system on which it was intended to be installed on.

To return a Technimount product,

- A Return of Merchandise Authorized (RMA) must be requested and approved by Technimount prior to returning the product.
- Products must be returned undamaged and in its original packaging, appropriately identified with the approved RMA number. Returns will not be approved on a modified or damaged item.
- Charges may apply if the package received is damaged or items are missing.
- Purchaser is responsible for a restocking fee (refer to Table 1 on page 8).

**Table 1:** Restocking fees

RESTOCKING FEES	
Prior to thirty (30) days	10%
Prior to forty-five (45) days	25%
Prior to sixty (60) days	30%

For any manufacturing defect, refer to the conditions within the warranty policy or contact Customer Service at [customerservice@technimount.com](mailto:customerservice@technimount.com) for additional information.

#### 1.4.3. Return of Material Authorization (RMA)

The Technimount Customer Service department is responsible for all merchandise returns and will provide a Return of Merchandise Authorization (RMA) number, upon approval. The RMA must be printed and placed on the returned merchandise. Technimount reserves the right to charge shipping and restocking fees (refer to Table 1) for the returned items. Special, modified, or discontinued items are not subject to returns.

#### 1.4.4. Claim Process

Upon reception of the returned merchandise, a thorough inspection will be performed. If the merchandise is compliant with the return policy or it is found that the product is defective, Technimount will take corrective actions and close the claim. If, however, it is found that the product is not defective, but rather misused or abused, the product will not be covered by the warranty. Details of our findings and conclusions will be provided shortly thereafter. To submit a claim, contact Customer Service at [customerservice@technimount.com](mailto:customerservice@technimount.com) to obtain a Return of Material Authorization (RMA) form and return instructions.



## 2. General Safety Guidelines

Always read and abide by all the safety guidelines identified within this document. Pictograms, safety symbols and labels are used to alert the user to a potentially hazardous situation, which, if not avoided, may result in minor or moderate injury to the patients or EMS and clinical personnel, or damage to the product. This includes the special care necessary for the safe and effective use of the Technimount product to avoid damage that may occur from use or misuse. The terms "Warning" and "Caution" herein carry special meaning and should be carefully reviewed.

**WARNING** – Indicates a hazardous situation that, if not avoided, could result in death or serious injury.

**CAUTION** – Indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.

### 2.1. Symbols and Definitions



**WARNING – Hand Crush/Pinch Point**

Indicates an area where mechanical components could move toward each other and might result in a potential crush/pinch hazard.



**WARNING – Do Not Step**

Indicates an area where there is potential risk of tipping if the user steps, stands, sits or rests his/her foot that could result in serious injuries to the patients or EMS and clinical personnel, or damage to the product.



**CAUTION – Safe Working Load (SWL)/Load Balance**

Indicates the total maximum charge for a safe use of the product.



**CAUTION – Safe Handling and Operation**

Alerts the reader to pay special attention to the recommendations for safe use of the product, and of potentially hazardous situations that could result in minor injuries to the patients or EMS and clinical personnel. This includes the special care necessary for the safe and effective use of the product to avoid damage that may occur from use or misuse.



**CAUTION – Safe Practice**

Alerts the reader to pay special attention to the recommendations and methods outlining how to safely operate the product to minimize risks to the patients, EMS and clinical personnel and the product.



**CAUTION – General Mandatory Action**

Call for action. Alerts the reader to potential risk to the patients or EMS and clinical personnel not following the mandatory action specified by the supplementary sign.



**CAUTION – Follow Instructions for Use**

Call for action. Reminds the reader to consult the user manual for information.




**CAUTION – Two (2) Person Lift**

Heavy load. Alerts the reader to a two (2) person lift carrying technique recommendation based on the weight and/or size of the product.

## 2.2. Labels

Labelling on the surface of the Technimount product, quickly identify potential risks and provides information to the user. A manufacturing label, including the serial number (Figure 1), can be seen on the Technimount product.



1. Manufacturing label

**Figure 1:** Location of the manufacturing label (bottom of the mounting system)

### 2.3. Safety Measures

Carefully read all the safety measures herein before operating the Technimount product, relay to EMS and clinical personnel during training, and include in your existing protocols.

More specific safety measures intended for biomedical technicians (or equivalent) relating to the safety checks and conditioned-based maintenance can be found in « Annex IV Maintenance » on page 27.



#### **WARNING – Risk of Injury**

- **Do not** use the mounting system if the mounting screws are loose or missing, to prevent undue risk to the medical device, patients, and EMS and clinical personnel.
- Always use compatible medical devices and accessories when applicable, to avoid unpredictable functioning resulting injury to the patients or EMS and clinical personnel. Refer to the « Technical Specifications » on page 12 for compatibilities.
- Improper use of the Technimount product may damage it or cause injury to the patients or EMS and clinical personnel.
- If any serious incident occurs with the mounting system, immediately stop using the product, report this incident to Technical Support at [technicalsupport@technimount.com](mailto:technicalsupport@technimount.com) and the applicable regulatory agency.



#### **CAUTION – Safe Practice**

- Practice safely operating the mounting system until the manipulations have been perfected, before use with patients. Improper use of a Technimount product may damage it or cause injury to the patients or EMS and clinical personnel.
- Regulations and standards for safety are the sole responsibility of the end user. Ensure that the installation specifications meet the local and regional compliance requirements before use.
- Refer to your protocols and the user documentation provided with each specific medical device for the safety guidelines and safe use.



#### **CAUTION – Safe Handling and Operation**

Always ensure that the Wall Mount Pro - FL is secured to the vertical surface before installing a mounting system and/or medical device in its Standard Surface Base, to avoid risks of damage, equipment falling, or injuries to the patients or EMS and clinical personnel.



#### **CAUTION – Safe Working Load (SWL)/Load Balance**

**Do not** overload the mounting system to avoid tipping incidents or risks of collapsing. For specifications, contact Customer Service at [customerservice@technimount.com](mailto:customerservice@technimount.com).



#### **CAUTION – Follow the Instruction for Use**


- Always read and abide by all the safety guidelines identified, as well as follow instructions provided within the user manual of the Technimount product.
- The vertical mounting system may contain optional medical equipment and accessories. Refer to their specific user documentation for the safety guidelines and safe use.

### 3. Technical Specifications

<b>Product Name</b>	Wall Mount Pro - FL
<b>Description</b>	Mounting system designed to aid trained EMS and clinical personnel to hold and secure on vertical surfaces in airplanes and helicopters, medical devices and accessories
<b>Product Code</b>	600-11-WMA-FL
<b>Operating Environment</b>	EMS/CCT (air)
<b>Compliance</b>	Designed to comply with FAA, 14 CFR § 23, 25, 27 & 29 (applicable sections)
<b>Expected Service Life</b>	5 years
<b>Compatible Stretcher</b>	N/A
<b>Compatible Mounting System</b>	<ul style="list-style-type: none"> <li>- Mounting Plate - FL</li> <li>- Airplane/Helicopter vertical surface</li> </ul>
<b>Compatible Medical Devices/ Accessories</b>	N/A
<b>Dimensions (W X D X H)</b>	7.5 in. X 10.09 in. X 7.19 in. (19.05 cm X 25.63 cm X 18.26 cm)
<b>Weight</b>	<ul style="list-style-type: none"> <li>- Wall Mount Pro - FL: 4.3 lb (1.95 kg), including the Standard Surface Base</li> <li>- Standard Surface Base: Refer to the user documentation</li> </ul>
<b>Composition</b>	<ul style="list-style-type: none"> <li>- Wall Mount Pro - FL: Aluminum</li> <li>- Standard Surface Base: Refer to the user documentation</li> </ul>
<b>Total Safe Working Load (SWL)</b>	For specifications, contact Customer Service at <a href="mailto:customerservice@technimount.com">customerservice@technimount.com</a>
<b>Operating Temperature</b>	- 31° F to 113° F (- 35° C to 45° C)
<b>Cleaning Solutions</b>	<ul style="list-style-type: none"> <li>- Oxivir®, 5% Hydrogen Peroxide with Peracetic Acid (AHP)</li> <li>- Lavo® 12, 10 000 ppm Sodium Hypochlorite</li> <li>- TNT-100, 5% Quaternary Ammonium Compound</li> <li>- Spectro-Sept, 5% Ethyl Alcohol</li> <li>- Spectrol, 5% EDTA salt</li> </ul>
<b>Options</b>	N/A


#### 4. Wall Mount Pro - FL Orientation Diagrams

**NOTE :** The orientations referenced herein are from the EMS and clinical personnel standpoint, when facing the mounting system on the vertical surface of the airplane or helicopter.




**Figure 2:** Wall Mount Pro - FL orientation diagram

5. Wall Mount Pro - FL Illustrated Parts



- |   |  |
|---|--|
| 1. Standard Surface Base                        | 3. Wall mounting installation screw holes (4X) |
| 2. Quick release button (Standard Surface Base) | 4. Wall Mount Pro - FL                         |

**Figure 3:** Wall Mount Pro - FL components



- |  |                                |
|--|--------------------------------|
| 1. Medical device mounted in a Bracket Pro Serie mounting system | 3. Disc of the mounting system |
| 2. Standard Surface Base   | 4. Wall Mount Pro - FL         |

**Figure 4:** Wall Mount Pro - FL with medical devices and accessories


## 6. Operate the Wall Mount Pro - FL

The contents in this section is intended for EMS and clinical personnel who are proficient, have received the required training and passed the skills assessment, therefore sufficiently knowledgeable to safely operate the mounting system.

**NOTE :** Illustrations throughout this user documentation show the Wall Mount Pro - FL and a medical device secured in a Bracket Pro Serie mounting system for comprehension purposes, but the same instructions apply for any compatible mounting systems and/or medical devices. Refer to the « Technical Specifications » on page 12 for the compatibilities and to the user documentation of each component for the safety guidelines and safe use.

### 6.1. Install a Mounting System on the Wall Mount Pro - FL


1. Align and insert the disc of the mounting system in the base of the Wall Mount Pro - FL horizontally (Figure 5 A).
2. Push the mounting system all the way back in the base until it is locked (Figure 5 B).



**Figure 5:** Installing the mounting system on the Wall Mount Pro - FL



3. Ensure that the mounting system is secured. If the disc stays in the base after the verification, it is properly secured.
4. Turn the mounting system clockwise or counterclockwise (Figure 6), to the desired position when needed.




**Figure 6:** Operating the mounting system on the Wall Mount Pro - FL

The installation of the mounting system on the Wall Mount Pro - FL is complete.

## 6.2. Remove a Mounting System from the Wall Mount Pro - FL

1. Press and hold the quick release button of the base (Figure 7 A), then pull the mounting system horizontally out of the base (Figure 7 B).



**Figure 7:** Removing the mounting system from the Wall Mount Pro - FL

2. Set aside the mounting system aside on a flat and clean surface or store as per your internal protocols.

The removal of the mounting system from the Wall Mount Pro - FL is complete.

## Annex I EMS and clinical personnel Skills Assessment

Following training, a skills assessment should be given to each member of the EMS and clinical personnel to ensure they have fully comprehended the labelling, warnings and cautions, potential risks, safe practices and proper operating procedures needed to safely and effectively use the mounting system. Consider adding the following to your internal training protocols.

Trainee name: \_\_\_\_\_ Unit: \_\_\_\_\_

Assessor name: \_\_\_\_\_ Date: \_\_\_\_\_

### EMS AND CLINICAL PERSONNEL SKILLS ASSESSMENT

SKILL CRITERIA	PASSED	FAILED
<b>Safety Measures</b>		
- Knows not to use the Wall Mount Pro - FL if the mounting screws are loose or missing.	<input type="checkbox"/>	<input type="checkbox"/>
- Knows to ensure that the medical device or accessory is secured in the mounting system before it is moved.	<input type="checkbox"/>	<input type="checkbox"/>
- Knows not to overload the mounting system and its components.	<input type="checkbox"/>	<input type="checkbox"/>
<b>Operation</b>		
- Able to install/remove a mounting system on/from the Wall Mount Pro - FL.	<input type="checkbox"/>	<input type="checkbox"/>
- Able to operate a mounting system on the Wall Mount Pro - FL.	<input type="checkbox"/>	<input type="checkbox"/>
- Has practiced safely operating the Wall Mount Pro - FL, has perfected the manipulations and has acquired the required skill level to safely use with patient.	<input type="checkbox"/>	<input type="checkbox"/>



## Annex II Unpack the Wall Mount Pro - FL

Unpacking should be reserved for biomedical technicians (or equivalent) who have extensive mechanical experience, and advanced skill level.

1. Inspect the shipping box(es) for signs of damage before accepting shipment. Take pictures and report them promptly when applicable.
2. Move the shipping box(es) to the location of the installation.
3. Open the shipping box(es).
4. Unpack the box(es) and ensure that all shipping and packaging materials have been properly removed, prior to installation.

**NOTE :** Keep all packaging material for future use.

5. Identify all the components and hardware included for the installation when applicable, then set aside. Refer to Annex III on page 23 for the required parts.
6. Inspect the items for signs of damage. Take pictures and report them promptly when applicable.



## Annex III Install the Wall Mount Pro - FL



### **WARNING** – General Warning

Always ensure that the vertical surface can withstand the weight of the Wall Mount Pro - FL with the medical device and accessories. Under no circumstances will Technimount be held liable hereunder for incidental or consequential damages, arising from a non-conforming installation.

The contents in this section is intended for biomedical technicians (or equivalent) who have extensive mechanical experience, and advanced skill level and have read the « Safety Measures » on page 11.

**NOTE :** The following installation procedure illustrates the Wall Mount Pro - FL being installed on a Mounting Plate - FL, but the same instructions apply for any compatible vertical surface.


### Required Installation Time

5 minutes

### Required Tools

- 1/4-28 MS24693, 100 deg. stainless steel screws (4X; User provided. The length will depend on the installation surface)
- 1/4 in.-28 ASTM F594C stainless steel hex nut or equivalent are recommended if bolting the Mounting Plate - FL on the vertical surface of the airplane/helicopter (4X; User provided)
- Torque wrench
- 7/16 wrench
- Screwdriver (Depends on the hardware)
- Medium strength thread lock adhesive (🔧)

Required Parts




- |  |                        |
|--|------------------------|
| 1. Standard Surface Base (factory installed)                           | 3. Wall Mount Pro - FL |
| 2. 1/4-28 MS24693, 100 deg. stainless steel screws (4X; user provided) |                        |

**Figure 8:** Wall Mount Pro - FL assembly parts

1. Ensure that the chosen installation surface meets the safety and installation requirements. Refer to the « Safety Measures » on page 11 and to the « Technical Specifications » on page 12 if needed.
2. For installation directly on the vertical surface of the airplane/helicopter, drill/tap holes on the installation surface if needed.
3. For installation on a Mounting Plate - FL or directly on a vertical surface, position the Wall Mount Pro - FL on the installation surface, then secure the mounting system using four (4) 1/4-28 MS24693, 100 deg. stainless steel screws coated with medium strength thread lock adhesive and a screwdriver (Figure 9).

**NOTE :** Apply a torque of 65-73 in.-lb to each screw.



**Figure 9:** Installing the Wall Mount Pro - FL on a Mounting Plate - FL



4. For installations in a bolted vertical surface, bolt the Wall Mount Pro - FL using four (4) 1/4-28 MS24693, 100 deg. stainless steel screws coated with medium strength thread lock adhesive (the length of the screws will depend on the installation surface), four (4) 1/4 in.-28 ASTM F594C stainless steel hex nut or equivalent and a screwdriver and wrench.

**NOTE :** Apply a torque of 65-73 in.-lb to each screw.

5. Ensure that the Wall Mount Pro - FL is properly secured on the installation surface. It is properly secured if it does not move.

The installation of the Wall Mount Pro - FL on the vertical surface of the airplane/helicopter or Mounting Plate - FL is complete.



## Annex IV Maintenance

Safety checks and condition-based maintenance should be carried out by biomedical technicians (or equivalent) who have extensive mechanical experience, and advanced skill level and have read all the « Safety Measures » on page 11, and the maintenance specific safety measures listed below.

Factors such as weather, environment, geographical location and individual usage will necessitate different needs. For the maintenance of the Wall Mount Pro - FL, follow the guidelines listed herein and in accordance with your service's current maintenance practices and protocols. Please contact Technical Support at [techsupport@technimount.com](mailto:techsupport@technimount.com) for replacement parts or repair related issues, if needed.



### **WARNING – General Warning**

- **Do not** perform safety checks or condition-based maintenance before having read the entire content of the user manual, gained in-depth knowledge and product comprehension, and familiarized yourself with the standards and guidelines.
- Safety checks and a condition-based maintenance plan are required and should be established for all Technimount products.
- Perform the safety checks and maintenance operations as described herein. Failing to follow the recommended maintenance plan or its guidelines could cause premature damage to the product.
- Use only Technimount parts, maintenance procedures, cleaning solutions and lubricants, as described herein. Using unapproved modified parts or procedures for the maintenance of the Technimount product may cause the system to be unstable and could cause injury to the patients or EMS and clinical personnel and void the product warranty.
- Replace damaged or worn-out parts if past their expected service life or when damaged (refer to « Annex V Replacement Parts/Kits » on page 33). Recycle damaged parts or dispose according to the environmental laws that apply to your jurisdiction and consult the Safety Data Sheets (SDS).



### **CAUTION – Safe Handling and Operation**

- **Do not** use powered tools to screw the hardware during installation, as there is a potential risk of damage to the threads.
- **Do not** steam clean or use ultrasonic cleaners on the system or any of its components.
- **Do not** immerse the metal parts/components in water.
- To spot clean, the maximum water temperature should not exceed 180° F/82° C. The maximum water pressure should not exceed 1500 psi/103.5 BAR. If using a high pressure washer, the pressure nozzle must be kept a minimum of 24 in. (61 cm) from the unit.
- When cleaning, always use appropriate Personal Protection Equipment (PPE) based on established protocols (e.g., gloves, eyewear, etc.).



### **CAUTION – Corrosion**

- Always rinse and dry the mounting system properly after using cleaning products. Certain types of cleaners may leave a corrosive residue on the surface of the product and could cause the premature corrosion of critical components. Refer to the product Safety Data Sheets (SDS) for chemical information or handling, storage and emergency measures in case of accident.
- Dispose of corrosive wastes according to the environmental laws that apply to your jurisdiction and consult the Safety Data Sheets (SDS).

**CAUTION – Follow Instructions for Use**

Always read and abide by all the safety guidelines identified, as well as follow instructions provided by the manufacturer of the cleaning product.

**Maintenance Frequency**

- Safety checks and the condition-based maintenance should be performed minimally every month or as frequently needed, to prolong the longevity of the mounting system in optimal conditions.
- Decontaminate the mounting system as recommended in your internal protocols, as well as the regulations and standards in virtue of the infection prevention and control procedures.

**Required Tools**

- Clean dry cloths
- Soft brush
- Pressure washer
- Cleaning solutions
- Medium strength thread lock adhesive (🔩)
- Phillips-head screwdriver #3
- Torque wrench

**Tested Cleaning Solutions**

- Oxivir, 5% Hydrogen Peroxide with Peracetic Acid (AHP)
- Lavo 12, 10 000 ppm Sodium Hypochlorite
- TNT-100, 5% Quaternary Ammonium Compound
- Spectro-Sept, 5% Ethyl Alcohol
- Spectrol, 5% EDTA salt

### Maintenance Plan

**NOTE :** In case of a non-conformity, stop using the product and contact Technical Support at [techsupport@technimount.com](mailto:techsupport@technimount.com) immediately for a remedial action plan.

**NOTE :** Always keep records of your maintenance activities and immediately remove defective or expired products from your inventory.

MAINTENANCE PLAN	COMPLIANT	
SAFETY CHECKS	YES	NO

#### Wall Mount Pro - FL (Figure 10)

- Visually inspect all the components of the mounting system to ensure there is no damage or chemical attack, that the hardware is in good condition and there are no loose screws:
  - Wall mounting installation screws
  - Standard Surface Base and quick release button. Refer to the user documentation.
  - If there is damage to the components, remove the product from circulation and contact Technical Support immediately for a remedial action plan.
  - If there are traces of chemical attack, follow the conditioned-based maintenance herein.
  - If the hardware is not in good condition, replace it. Contact Technical Support if needed.
  - If the hardware is loose, apply medium strength thread lock adhesive and tighten using a Phillips screwdriver.

CONDITION-BASED MAINTENANCE	YES	NO
-----------------------------	-----	----

Following the safety checks,

#### Clean the Wall Mount Pro - FL

1. Remove the excess dirt using a clean cloth, if needed.
2. Remove the contaminants using a pressure washer or as recommended in your internal protocols and control procedures.
3. Clean using a cloth and cleaning solution.
4. Spot clean stains by applying the solution directly on the stain and let sit on the surface, if needed.

**NOTE :** Avoid over saturation and ensure that the product does not sit on the surface of the mounting system longer than recommended by the cleaner's manufacturer.


5. Thoroughly rinse the solution with a clean cloth dampened with lukewarm water, then dry all the components using a clean cloth before returning to service.

Comments and observations following the Safety Checks and Condition-Based Maintenance:

Maintenance plan completed on (dd/mm/yyyy):

Maintenance plan completed by:

Illustrated Inspection Points



- |   |  |
|---|--|
| 1. Standard Surface Base                        | 3. Wall mounting installation screw holes (4X) |
| 2. Quick release button (Standard Surface Base) | 4. Wall Mount Pro - FL                         |

**Figure 10:** Illustrated inspection points





## Annex V Replacement Parts/Kits

Technimount reserves the right to change part numbers and products without notice. Please contact Customer Service at [customerservice@technimount.com](mailto:customerservice@technimount.com) to ensure product options and availability, or Technical Support at [techsupport@technimount.com](mailto:techsupport@technimount.com) for replacement parts/kits or repair related issues.

PART/KIT NUMBER	PART/KIT DESCRIPTION
N/A	N/A



**TECHNIMOUNT**  
EMS<sup>®</sup>

SAFETY AND FLEXIBILITY  
WHERE IT MATTERS MOST